


SURGICAL GOWN

FDA 510(k) CLEARED

Choose SAFETY for your gown.



MEET THE REQUIREMENTS OF EN & ANSI/AAMI STANDARD FOR SURGICAL GOWN

EN13795:2011 Surgical drapes, gowns, and clean air suits, used as medical devices for patients, clinical staff, and equipment — General requirements for manufacturers, processors and products, test methods, performance requirements and performance levels.

ANSI/AAMI PB70 : 2012 Liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities.

FEATURES & BENEFITS

SAFETY CONCEPTS

- Double neck closure (Velcro and ties)
- Soft stretchy knit cuff is long enough to avoid contamination
- Multi folded hand towel intended for a sterile field
- Waist strap (front tie) attached within the critical zone

CONVENIENT DETAILS

- Specification stamped by case type for efficient identification
- Minimal touch to open the gown
- Easy to doff and discard (moderate strength of the gown ties)

FUNCTIONAL FABRIC (APPROPRIATE TO INTERNATIONAL STANDARDS)

- Meet the AATCC42 and AATCC127 AAMI Level 3 requirements within the critical zones.
- Meet EU Standard EN13795 (bacteria infiltration, lint test, particulate test, etc.)

GOWN SEAMS


- Waist straps by ultrasonic adhesion
- Two Types of Sleeve Seams: Sewing/Ultrasonic
- Protect from fluid exposure

HIGH ABSORBENT HAND TOWEL

- Economical, highly absorbent and latex-free
- Low lint reduces particulates

SMMMS Structure of Fabric

SMMMS(5Layers) Structure of Fabric



Sejong Healthcare’s premier surgical gown, made of high quality spun-melt nonwoven fabrics, offers the safety of any infection and pleasant wearing sensations for doctors and patients. Boasting excellent water repellent, blood resistant and antimicrobial properties, specially treated multi-layered nonwovens provides users with outstanding protection against leakage of blood, alcohol and any liquids and good air-permeability.

- **Excellent Protection**  
High Level of hydrostatic pressure, anti-alcohol, blood repellency and anti-bacteria
- **Breathable and Comfortable**
- **Exceptional strength and durability**
- **Low lint and dust**

High Water Repellency	Fabric of low lint and dust	High Resistance to Microbial Penetration	Anti-Statics	Anti-Alcohol
-----------------------	-----------------------------	--	--------------	--------------